

THE APPLICATION OF
RenuLife Violet Ray

High Frequency for

ASTHMA
BRONCHITIS
CATARRH (NOSE
AND THROAT)
COLDS
HAY FEVER

DEAFNESS AND EAR
DISEASES
SORE THROAT AND
THROAT DISEASES
WHOOPIING COUGH
PULMONARY AILMENTS

Asthma is a spasmodic contraction of wind pipe or bronchial tubes causing difficulty in breathing and a wheezing when breath is expelled. Authorities seem to agree that the primal cause is nervousness. There is evidence of inflammation in lining of wind pipe and smaller branches of bronchial tubes.

Both the inflamed local conditions and the disturbed nervous center causing the trouble should respond to High Frequency Violet Ray Treatments. The exterior of throat and chest is given strong current treatments. Ozone Generator is specified as well as surface treatments along spinal nerve centers.

Bronchitis—The inflammation that is centered in the bronchial tubes is called bronchitis. The area of inflammation, however, may extend to the windpipe and even to the lungs. Soreness in the top of the chest, coughing, fever, headache, are symptoms of the disease. A continuation of thorough treatment over surface of chest, throat and back with ozone inhalations is very effective.

Catarrh—When the mucous membrane, which lines all organs that open outside the body, becomes inflamed it is a condition known under the general name catarrh. This inflammation in the air passages, nasal and throat, is generally spoken of as catarrh. Violet Ray High Frequency unquestionably offers the catarrh sufferer more relief and benefit than can be expected from drugs. Both the exterior of the nose, treated with general applicator, the interior with special internal nose applicator, and inhalation from Ozone Generator for the triple attack on this dangerous disease.

Colds are an acute attack of catarrh. Treatment for colds is similar to that for catarrh and the results are usually rapid and satisfactory.

Hay Fever is a persistent irritation of the mucous lining of the nostrils and throat which recur in a person at a certain period of each year. Sneezing, watering nose and eyes and smarting throat and nose are some of the symptoms. Many reports have been made of the favorable influence of High Frequency in these cases. Direct application to the nose with the nasal electrode and mild current, also treatment of the nose externally. Inhalations of Ozone are of even greater importance.

Deafness and Ear Diseases—Deafness may be due to any one of several causes and the manner of treating for best results depends on the cause of the trouble. The ear drum may be diseased or it may have been punctured. It may be due to the disease of the nerves or small bones in ear. Inflammation of the eustachian tube which leads from throat to middle ear will cause deafness. For deafness and ear trouble High Frequency is strongly recommended. The special electrode for inserting in the ear should be used. A sensation of warmth is produced and congestion, inflammation and pain are almost immediately relieved.

Sore Throat and Throat Diseases—The application of High Frequency Violet Ray to the exterior of the throat is one of the most effective means for treating such local troubles as tonsillitis, quinsy and sore throat. High Frequency brings a fresh supply of blood to the affected parts, relieves congestion and pain.

• **Whooping Cough**—The efficacy of Ozone generated by High Frequency and taken through the Ozone Generator in treating whooping cough has been attested by many authorities. Says L'Abbe: "Ozone is the remedy par excellence for whooping cough."

Pulmonary Ailments—Ozone is especially effective in consumption. This is due to the fact that Ozone is twice as powerful in oxidizing properties as ordinary air. For all diseases of the air passages the use of Ozone is strongly recommended.

OZONE

Ozone is simply a concentrated form of Oxygen or pure air. Without oxygen we cannot live. When we breathe we fill the lungs full of pure air or oxygen. This oxygen passes through the walls of the lungs at the same time passing through and purifying the blood which is flowing through the walls of the lungs and which comes there from every part of the body to be purified.

It is the lack of pure air to purify our blood that causes us to be overcome with disease. Being cooped up in an office or factory all day long is responsible for most of our diseases and particularly diseases of the respiratory organs such as Catarrh, Bronchitis, Pneumonia, Asthma, Tuberculosis and Colds.

To overcome these diseases most doctors advise a trip to the mountains, the seashore or the pine forests in order to obtain the desired ozone or oxygen. The **RENULIFE OZONE GENERATOR** generates ozone or oxygen in sufficient quantities to prevent all this expense if used in your own home.

The **RENULIFE OZONE GENERATOR** is an attachment which fits in the handle of any of our **VIOLET RAY** Generators and delivers Ozone in large quantities for inhalation purposes. The Ozone penetrates every cell in the lungs, purifying the blood, destroying disease germs and soothing the inflamed tissues. For Asthma, Catarrh, Hay Fever, Coughs, Colds, Bronchitis, Tuberculosis, etc., there is no better treatment known. As one of our most noted surgeons puts it, "Ozone is essential wherever pure air is essential and should be used wherever any disease is present."



Ozone Generator No. 25 \$9.50

INSTRUCTIONS FOR OPERATING

Ozone generators are packed with nose piece attached; included separately are a bottle of pine needle oil and a glass funnel.

First, see that there is a small quantity (about a teaspoonful) of pine needle oil in upper chamber of ozone generator. To put the oil in chamber, use funnel. Detach nose piece from rubber coupling, and insert funnel in its place.

A teaspoonful of oil will do for a considerable period.

Insert stem of ozone generator in handle the same as any other electrode. **Grasp holder of ozone generator with left hand**, as shown in illustration; turn on medium current, and agitate bulb with right hand. Ozone filtered through pine needle oil is produced. A pleasing, purifying, powerful agency.

THE APPLICATION OF
RenuLife Violet Ray

High Frequency for

ABSCESSSES	ECZEMA AND SKIN
BALDNESS, FALLING	DISEASES
HAIR (ALOPECIA)	PIMPLES (ACNE)
DANDRUFF	SCARS
BLACKHEADS	WARTS
BRUISES	WRINKLES
BUNIONS	

Abscesses—High Frequency Violet Rays are specified for treatment of abscesses, as purifying and stimulating effect destroys the diseased condition.

Baldness, Falling Hair (Alopecia), Dandruff—High Frequency is an ideal remedy for loss of hair, which results from faulty nutrition. Two daily treatments of five minutes' duration taken regularly and persistently almost invariably give satisfactory results.

Dandruff, gray hair, itching scalp are likewise successfully treated. Violet Rays are the best hair tonic procurable, as they actually revitalize the hair roots. In many cases gray hair has been restored to its original color.

Blackheads—This annoying skin disease is caused by faulty functioning of the skin. Violet Ray treatments not only will remove the blackheads, but, by stimulating the skin and bringing blood to the treated area, establish a healthful condition, preventing a recurrence of the trouble.

Bruises—Applied to the bruised or swollen area, the Violet Ray generally brings immediate results. Best results are accomplished by using high frequency lubricant, which is driven into the tissue by the application of the Violet Ray.

Bunions—This painful affliction of the big toe joint is greatly relieved by Violet Ray treatments. The inflammation subsides and soreness disappears.

Eczema and Skin Diseases—The reaction to High Frequency Violet Ray on the skin is such

that diseases, both of germ origin and from malnutrition, are successfully combatted. Reports of cure by High Frequency currents of eczema are numerous.

Pimples (Acne)—When the skin is diseased, High Frequency Violet Ray, applied through vacuum in a shower of sparks over the affected area, is very effective.

Scars—Vacuum tube sparks have a softening effect on scar tissue. For some scars caustic fulguration spark should be used which will destroy the scar tissue and pliable new tissue will take its place.

Warts and Moles—Beauty parlors regularly employ the Violet Ray, applied through a fulguration electrode, to remove warts, moles and other facial blemishes. The spot is treated with a hot spark which sears it. A crust or scab forms after several treatments which later separates, leaving no scar. Sometimes it requires several treatments to entirely eradicate the growth.

Wrinkles—These are generally caused either by flabby facial muscles or by the habitual use of one set of facial muscles, causing lines to appear on the face. High Frequency Violet Ray revitalizes muscles and smooths out the wrinkles. Every beauty parlor makes use of this treatment; it can be just as effectively applied at home and be taken more often.

THE APPLICATION OF
RenuLife Violet Ray

High Frequency for

ANEMIA	FEMALE COMPLAINTS
ARTERIOSCLEROSIS	GOUT
(HARDENING OF	HEADACHE
ARTERIES)	NERVOUS AFFECTIONS
CIRCULATORY	PILES
DISORDERS	PROSTATIC DISEASES
CONSTIPATION	

Anemia—High Frequency Violet Ray currents increase the oxygenation of the blood. Their employment in anemia, therefore, is quite effective. They may be used without interference with iron tonics, if desired.

Ozone inhalations are of the utmost value. This greatly facilitates absorption of oxygen by the blood resulting in a rapid increase of hemoglobin.

Arteriosclerosis (Hardening Arteries) (High Blood Pressure)—The use of Violet Ray High Frequency home treatments for these ailments is beneficial, but the sufferer should place himself in the hands of an electro-therapist, as a general systemic thorough treatment with diet instructions are required by the seriousness of the affliction. The use of High Frequency will lower the blood pressure.

Circulatory Disorders—No surer means of stimulating and accelerating circulation is known than High Frequency Violet Ray treatments. At the same time the blood itself is purified. Deposits and impediments dissolved. Cold hands or feet indicate the absence of rich, warm blood. This chilliness is quickly overcome by applications of Violet Ray to the parts affected. At the same time the blood is built up and eventually normal circulation established.

Constipation—High frequency may be used advantageously for constipation, application being made over the abdomen. Also, special rectal electrode for internal treatment.

Female Complaints—Very gratifying results secured in such cases. Full instructions for treatment are given in our treatment chart.

Gout—Gout is usually the result of excessive eating of rich food, or drinking of alcoholic beverages. Piercing pains, causing great agony, are suffered in the joint of the great toe. The diet condition should be immediately attended to and local treatments of High Frequency over the area of pain will give quick and sure results.

Headache—High Frequency Violet Ray treatment is very effective for headache. Only temporary relief may be expected where it is the reflex of serious organic diseases. These causes should be determined and also treated. Not only the head itself, but the spine, neck cord, and stomach should be treated with High Frequency, depending upon the cause, of which a headache is the symptom.

Nervous Affections—Numerous cases of nervous exhaustion have been reported cured by High Frequency. The general systemic treatment through the metal electrode, held in the hand, is quieting and relieves nerve tension.

Ozone inhaled is a very desirable adjunct. The great nerves of spine and neck should also be treated.

Piles—The insulated electrode should be used. High Frequency current has proved most efficacious, whether used alone or in connection with some other remedy.

Prostatic Diseases—The results following the high frequency treatment of all forms of prostatic diseases are extraordinary. Authorities report complete cures in a very large percentage of cases. The vacuum tube, especially designed for such treatments, should be used.

THE APPLICATION OF
RenuLife Violet Ray

High Frequency for

LAMENESS (LOCO- MOTOR ATAXIA)	NEURITIS
INFANTILE PARALYSIS	PARALYSIS
LUMBAGO	RHEUMATISM
NEURALGIA	SCIATICA

Infantile Paralysis—Violet Rays are so mild and gentle in action that the smallest infant may be treated without harmful results. Infantile paralysis caused by accident or disease is a natural condition for Violet Ray treatment. The spine, as well as muscles affected, should be treated. The opportunity to give daily treatments with our instrument makes its use especially advantageous.

Lameness (Locomotor Ataxia)—Relief and benefit follows the use of high frequency in such cases.

Lumbago—Immediate relief follows the use of High Frequency for lumbago, which is a form of muscular rheumatism. The application should be made to the affected area with the patient changing his position often. Treatments may be taken as often as necessary and of ten minutes' duration until relief is given.

Neuralgia—The old definition of neuralgia was "the cry of a starved nerve for blood." This essentially indicates the use of High Frequency Violet Ray. No other agency is so useful in stimulating the circulation or bringing to the area increased blood supply. At the same time, it purifies the blood.

Neuritis—Neuritis is the inflamed condition of the nerve, causing pain. Any nerve may be involved, but frequently the nerve of shoulder and arm is affected. High Frequency is a positive relief and almost invariably a cure for this condition. Treatment should be in direct contact with the skin. At first the pain may be aggra-

vated, but the results are so satisfactory that the slight initial unpleasantness may well be borne in view of the great benefits.

Paralysis—High Frequency is a most important treatment for paralysis. Often its persistent use will restore, partially or completely, the organ, limb, or set of muscles. Regular treatments over a considerable period of time are necessary in more stubborn cases.

Rheumatism—High Frequency Violet Rays are of exceptional value in muscular and in chronic and articular rheumatism. Pain is quickly relieved. Sometimes it is necessary to continue the treatment of ten minutes' duration several times a day. Generally the treatments should be through clothing or several layers of towel in order to get the benefit of the warm sparks. Not only is the pain relieved, but high frequency furnishes elimination and nutrition demanded in such cases. The use of the Ozone Generator will be found of great value in addition to the body applications.

Sciatica—This is an affection of the nerve trunk line which runs down the thigh from hip to knee. They lie deep under the muscle and the penetrative force of the High Frequency Violet Ray reaches the irritated area and furnishes a relief impossible to secure otherwise.

THE APPLICATION OF
RenuLife Violet Ray
High Frequency for

BRAIN FAG

GOITRE

INSOMNIA

PAIN IN ABDOMEN
AND CHEST

OBESITY

PYORRHEA

SPRAINS

TOOTHACHE

WEAK EYES

Brain Fag—Brain fag is a condition just one step removed from nervous prostration. High Frequency current and Ozone accomplish remarkable results. Treatments over spine and abdominal region and back of neck have a quieting, exhilarating effect most welcome to the sufferer.

Goitre—Gratifying results are secured in many cases by the use of High Frequency Violet Ray in treating goitre. The external throat electrode is so shaped that treatment of the goitre itself can be accomplished most effectively. Daily treatments persistently kept up are required.

Insomnia—There are few cases of insomnia that cannot be cured by High Frequency currents. The back of the head, neck, eyebrows, and eyes should be treated.

Obesity—Slow but satisfactory results are secured by High Frequency Violet Ray in reducing surplus fat. The diet must be looked after and daily treatments taken.

Pain In Abdomen and Chest—High Frequency current has a remarkable pain-relieving power. Application should be made at the seat of the pain until relief comes. Actual area treated glows with the warmth and cause of irritation is removed.

Pyorrhea—Many people troubled with pyorrhea are anemic. This would indicate a necessity of treating the blood as well as a local treatment of the diseased gums. Many dentists employ High Frequency Violet Ray in their work. Not only should the gums be treated with the special dental electrodes, but we strongly recommend

that the patient take inhalations from Ozone Generator. High Frequency Violet Ray is an antiseptic as well as stimulating.

Sprains—High Frequency Violet Ray removes soreness, stiffness from sprained muscles, tends to dispel the congestion, and is the pleasant and sure relief.

Toothache—In many cases of toothache the application of Violet Ray will give immediate relief. The cause, however, should be removed. This generally requires a dentist.

Weak Eyes—Special electrodes fitting the eye lids are best used in treating the eyes. A pleasant warmth is experienced, blood is brought to the affected area and soreness disappears. Muscular deterioration is stopped and gradual return to normal condition usually results.

RENULIFE ELECTRIC CO., INC.,
~~6th Floor Marquette Bldg.,~~
DETROIT, MICH.

*1304 Harper Str.,
Detroit, Mich.*

IMPORTANT

You Must Always Hold Handle at End Nearest Cord

Holding handle near electrode end may cause a discharge of current through handle, thus detracting from current delivered for treatment through electrode.



TO INSERT OR REMOVE ELECTRODE FROM HANDLE

1st—Grasp electrode with thumb and first finger of left hand and applicator handle in right hand as shown in upper photo.

2nd—Insert electrode into handle until electrode strikes electrode holder. (Slight resistance is felt).

3rd—Press electrode firmly into place.

TO REMOVE REVERSE THE ABOVE OPERATION

DON'T TWIST ELECTRODE while inserting or removing from handle

Renulife Electric Company

Detroit,

Michigan